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## DISPOSITIF MEDICAL

DELTA MED SpA Viadana (MN)	<b><u>MEDICAL DEVICES — TECHNICAL FORM</u></b>	Date:01/09/14 rev.10 Issued by: QA/QM Page 2 of 10
Family	Single and dual Entry I.V. Catheters	
Device	Invasive, surgical type medical device, non active, for short-term use (more than 60 minutes, less than 30 days) (REF. 93/42/EEC/MDD AND 47/2007 EEC/MDD – APPENDIX IX) — Disposable	
Classification	<b>Class II a</b> (REF. MDD 93/42/CEE – APPENDIX IX) – <b>Rule 7</b> <b>Class II</b> (RIF. CMDR SOR 98-282 – SCHEDULE 1) – <b>Rule 1</b> <i>UMDNS: 10727 – Catheters, Intravenous, Peripheral</i> <i>CND (D.M. 20/02/2007): C0101010102 – Aghi cannula senza valvola d'iniezione</i> <i>C0101010101-Aghi cannula con valvola di iniezione</i> <i>COD. CATEG. REF. prEN1874: 10 – Disposable devices</i>	
Definition	<b>Short-term peripheral venous access.</b> Device for peripheral intravenous catheterisation. Only in connection with other medical devices such as infusion set and/or needleless syringes containing a solution it permits the administration.	
Duration of contact	In order to reduce the risk of thrombophlebitis and bacterial invasion in the CVP <i>in situ</i> in the short term, it is recommended that the device stays in the vein 72-96 hours maximum. (Ref. Atlanta CDC)	

### LIST OF DISPOSABLE MEDICAL DEVICES, SUBJECT OF THIS TECHNICAL SHEET

**PRODUCT:** I.V. peripheral catheter (cannula needle), consisting of stainless steel needle, catheter in Polyurethane, with or without Latex-free injection valve. The peripheral venous catheter is commonly used in hospital and outside hospital environments to allow venous access in quick times.

DEVICE	CODE / TYPE
Single Entry I.V. Catheter in PUR (without wings) Catheter in Polyurethane(PUR) without Latex-free injection valve The various codes countermark the diameter of the single entry I.V. catheter, without Latex-free valve, in PUR.	3179022 / NEO DELTA VEN T 14G 32 MM
	3178522 / NEO DELTA VEN T 14G 45 MM
	3179322 / NEO DELTA VEN T 14G 50 MM
	3179622 / NEO DELTA VEN T 14G 51 MM
	3179722 / NEO DELTA VEN T 14G 70 MM
	3169522 / NEO DELTA VEN T 16G 30 MM
	3169022 / NEO DELTA VEN T 16G 32 MM
	3168522 / NEO DELTA VEN T 16G 45 MM
	3169422 / NEO DELTA VEN T 16G 50 MM
	3169822 / NEO DELTA VEN T 16G 51 MM
	3158522 / NEO DELTA VEN T 17G 45 MM
	3159822 / NEO DELTA VEN T 17G 51 MM
	3149322 / NEO DELTA VEN T 18G 30 MM
	3188522 / NEO DELTA VEN T 18G 32 MM
	3148522 / NEO DELTA VEN T 18G 45 MM
	3149022 / NEO DELTA VEN T 18G 48 MM
	3149822 / NEO DELTA VEN T 18G 51 MM
	3139522 / NEO DELTA VEN T 20G 25 MM
	3139222 / NEO DELTA VEN T 20G 30 MM
	3138522 / NEO DELTA VEN T 20G 32 MM
	3139622 / NEO DELTA VEN T 20G 45 MM
	3139322 / NEO DELTA VEN T 20G 48 MM
	3139422 / NEO DELTA VEN T 20G 50MM
3128522 / NEO DELTA VEN T 22G 25 MM	
3118522 / NEO DELTA VEN T 24G 19 MM	
3108522 / NEO DELTA VEN T 26G 19 MM	

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## DISPOSITIF MEDICAL

DELTA MED SpA Viadana (MN)	<b><u>MEDICAL DEVICES — TECHNICAL FORM</u></b>	Date:01/09/14 rev.10
		Issued by: QA/QM
		Page 3 of 10
Family	Single and dual Entry I.V. Catheters	
Device	Invasive, surgical type medical device, non active, for short-term use (more than 60 minutes, less than 30 days) (REF. 93/42/EEC/MDD AND 47/2007 EEC/MDD – APPENDIX IX) — Disposable	
Classification	<p>Class II a (REF. MDD 93/42/CEE – APPENDIX IX) – Rule 7</p> <p>Class II (RIF. CMDR SOR 98-282 – SCHEDULE 1) – Rule 1</p> <p>UMDNS: 10727 – Catheters, Intravenous, Peripheral</p> <p>CND (D.M. 20/02/2007): C0101010102 – Aghi cannula senza valvola d'iniezione</p> <p style="padding-left: 40px;">C0101010101-Aghi cannula con valvola di iniezione</p> <p>COD. CATEG. REF. prEN1874: 10 – Disposable devices</p>	
Definition	Short-term peripheral venous access. Device for peripheral intravenous catheterisation. Only in connection with other medical devices such as infusion set and/or needleless syringes containing a solution it permits the administration.	
Duration of contact	In order to reduce the risk of thrombophlebitis and bacterial invasion in the CVP <i>in situ</i> in the short term, it is recommended that the device stays in the vein 72-96 hours maximum. (Ref. Atlanta CDC)	

PHYSICAL CHARACTERISTICS	
CHARACTERISTIC	<p>Short-term peripheral venous access. Device for peripheral intravenous catheterisation (PIC). Only in connection with other medical devices such as infusion set and/or needleless syringes containing a solution it permits the administration.</p> <p>MOD.: NEO DELTA VEN.</p>
PERIPHERAL INTRAVENOUS INTRAVASCULAR CATHETER	Rx-opaque flexible tube with cannula, intended to be placed in peripheral venous access (distal end of the Catheter).It allows the channelling of infusion liquids during administration (continuous or intermittent) or the channelling of blood.
RX-OPAQUE LINE	<p>It allows the measurement of the incannulated catheter by X-rays.</p> <p>The radiopaque line is incorporated into the material (PUR) so that it is not possible the dispersion of BaSO<sub>4</sub> into patient blood.</p>
CATHETER HUB	<p>It allows the cannula to be extracted from the venous access, supporting it.</p> <p>It allows the catheter to be fixed and correctly maintain its position, through a female Luer-Lock connection, and to connect it to an infusion set for the continuous or intermittent administration of liquids and medicines or to connect it to a blood sampling syringe.</p> <p>In accordance with EN ISO 10555-5, the colour countermarks the diameter of the Gauge.</p>
INJECTION VALVE	<p>Variant present only in dual entry catheters.</p> <p>It allows the rapid administration of medicines without perforation.</p>
PORT CLOSURE	<p>Variant present only in dual entry catheters.</p> <p>It allows the closure of the secondary connection.</p>
INTRODUCER NEEDLE	<p>Catheter needle conforms to the EN ISO 10555-5 standard.</p> <p>Allows the catheter to be inserted into the peripheral venous access (venipuncture). Introducer needle with triple sharpening, atraumatic</p>
NEEDLE HUB	Allows the cannula needle to be extracted, supporting it.
VENT FITTING BODY WITH VENT FITTING FILTER	Transparent and equipped with a water-resistant stopper, allows both the rapid display of the venous return and the blocking of possible external blood contamination.
PROTECTOR	It allows the cannula needle to be protected.
MALE LUER LOCK CAP	It allows the momentary suspension of the treatment through the sealed closure of the nearest end (access) of the catheter hub.

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## DISPOSITIF MEDICAL

DELTA MED SpA Viadana (MN)	<b><i>MEDICAL DEVICES — TECHNICAL FORM</i></b>	Date:01/09/14 rev.10
		Issued by: QA/QM
		Page 4 of 10

Family	Single and dual Entry I.V. Catheters
Device	Invasive, surgical type medical device, non active, for short-term use (more than 60 minutes, less than 30 days) (REF. 93/42/EEC/MDD AND 47/2007 EEC/MDD – APPENDIX IX) — Disposable
Classification	Class II a (REF. MDD 93/42/CEE – APPENDIX IX) – Rule 7 Class II (RIF. CMDR SOR 98-282 – SCHEDULE 1) – Rule 1 UMDNS: 10727 – Catheters, Intravenous, Peripheral CND (D.M. 20/02/2007): C0101010102 – Aghi cannula senza valvola d'iniezione C0101010101-Aghi cannula con valvola di iniezione COD. CATEG. REF. prEN1874: 10 – Disposable devices
Definition	Short-term peripheral venous access. Device for peripheral intravenous catheterisation. Only in connection with other medical devices such as infusion set and/or needleless syringes containing a solution it permits the administration.
Duration of contact	In order to reduce the risk of thrombophlebitis and bacterial invasion in the CVP <i>in situ</i> in the short term, it is recommended that the device stays in the vein 72-96 hours maximum. (Ref. Atlanta CDC)

DESCRIPTION	REF. EN ISO 10555-5		CATHETER LENGTH (MM)	Ø EXT. CATHETER (MM)	Ø INT. CATHETER (MM)	FLOW (ML/MIN)	CATHETER MINIMUM BREAKING LOAD	Ø EXT. NEEDLE (MM)
	COLOUR CODE	GAUGE						
Dual Entry Catheter 26G	Violet	26	19 ±1	0.64	0.45	19	5N	0.4
Dual Entry Catheter 24G	Yellow	24	19 ±1	0.74	0.55	29	6N	0.5
Dual Entry Catheter 22G	Light blue	22	25 ±1	0.90	0.65	42	7N	0.6
Dual Entry Catheter 20G	Pink	20	32 ±1	1.00	0.75	59	8N	0.7
Dual Entry Catheter 18G	Green	18	32 ±1	1.30	0.95	103	11N	0.9
			45 ±1			96		
Dual Entry Catheter 17G	White	17	45 ±1	1.50	1.15	155	13N	1.1
Dual Entry Catheter 16G	Grey	16	45 ±1	1.75	1.35	225	16N	1.3
Dual Entry Catheter 14G	Orange	14	45 ±1	2.00	1.55	290	18N	1.5
Single Entry Catheter 26G	Violet	26	19 ±1	0.64	0.45	19	5N	0.4
Single Entry Catheter 24G	Yellow	24	19 ±1	0.74	0.55	29	6N	0.5
Single Entry Catheter 22G	Light blue	22	25 ±1	0.90	0.65	42	7N	0.6
Single Entry Catheter 20G	Pink	20	25 ±1	1.00	0.75	62	8N	0.7
			30 ±1			60		
			32 ±1			59		
			45 ±1			57		
			48 ±1			55		
			50 ±1			49		
Single Entry Catheter 18G	Green	18	30 ±1	1.30	0.95	110	11N	0.9
			32 ±1			103		
			45 ±1			96		
			48 ±1			94		
			51 ±1			93		
Single Entry Catheter 17G	White	17	45 ±1	1.50	1.15	155	13N	1.1
			51 ±1			152		
Single Entry Catheter 16G	Grey	16	30 ±1	1.75	1.35	233	16N	1.3
			32 ±1			230		
			45 ±1			225		
			50 ±1			223		
			51 ±1			222		
Single Entry Catheter 14G	Orange	14	32 ±1	2.00	1.55	300	18N	1.5
			45 ±1			290		
			50 ±1			288		
			51 ±1			287		
			70 ±1			280		

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## DISPOSITIF MEDICAL

DELTA MED SpA Viadana (MN)	<b><i>MEDICAL DEVICES — TECHNICAL FORM</i></b>	Date:01/09/14 rev.10 Issued by: QA/QM Page 5 of 10
Family	Single and dual Entry I.V. Catheters	
Device	Invasive, surgical type medical device, non active, for short-term use (more than 60 minutes, less than 30 days) (REF. 93/42/EEC/MDD AND 47/2007 EEC/MDD – APPENDIX IX) — Disposable	
Classification	<b>Class II a</b> (REF. MDD 93/42/CEE – APPENDIX IX) – <b>Rule 7</b> <b>Class II</b> (RIF. CMDR SOR 98-282 – SCHEDULE 1) – <b>Rule 1</b> <i>UMDNS: 10727 – Catheters, Intravenous, Peripheral</i> <i>CND (D.M. 20/02/2007): C0101010102 – Aghi cannula senza valvola d'iniezione</i> <i>C0101010101-Aghi cannula con valvola di iniezione</i> <i>COD. CATEG. REF. prEN1874: 10 – Disposable devices</i>	
Definition	<b>Short-term peripheral venous access.</b> Device for peripheral intravenous catheterisation. Only in connection with other medical devices such as infusion set and/or needleless syringes containing a solution it permits the administration.	
Duration of contact	In order to reduce the risk of thrombophlebitis and bacterial invasion in the CVP <i>in situ</i> in the short term, it is recommended that the device stays in the vein 72-96 hours maximum. (Ref. Atlanta CDC)	

THIS DATA ARE APPLICABLE ALSO TO ONE WAY CATHETER WITHOUT WINGS

### PRODUCTION METHOD

All the products are manufactured in compliance with the Legislative Decree 81/2008 and later updates as regards health and safety in the workplace and are manufactured in conformity with the company certified Quality System.

The parts used to manufacture the devices are produced using:

- moulding and extrusion (plastic parts of the catheters);
- drawing, grinding, polishing (needles)

They are checked in terms of physical (visual and dimensional) and functional aspect according to European and international standards and internal Delta Med procedures.

These devices are packed into the Clean Room. During production physical and functional checks is done in accordance with European and international standards and internal procedures. The finished products are then packaged inside the Clean Room in strips made from thermoformable film and thermosealed medical grade paper and then EO sterilised by the external sterilizing company.

*The products listed in this technical sheet are supplied in sterile form.*

### CLEAN ROOM ENVIRONMENT

The production is done in a Clean Room. The Clean Room, which is classified in accordance with European and international legislation EN ISO 14644-1 ISO CLASS 8, is periodically monitored, in "operational" state, with regards to the particulate and microbiological aspect and has the following properties:

- *Particulate contamination:* for particles  $\geq 0.5 \mu$   $\leq 3,520,000 \text{ m}^3$
- *Particulate contamination:* for particles  $\geq 5 \mu$   $\leq 29,300 \text{ m}^3$
- *Aerobiocontamination:*  $\leq 200 \text{ UFC/m}^3$  in accordance with Eu GMPs – Annex I (grade D)
- *Biocontamination of surfaces:*  $\leq 21000 \text{ UFC/m}^2$  (50 UFC for plate  $\varnothing 55 \text{ mm}$ ) in accordance with Eu GMPs – Annex I (grade D).

The access to the production area of materials and staff is done through a set of separate outer "lobbies". The hygienic rules and codes of conduct are observed by production staff and are shown in the access "lobbies" to the Clean Room.

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## DISPOSITIF MEDICAL

DELTA MED SpA Viadana (MN)	<b><u>MEDICAL DEVICES — TECHNICAL FORM</u></b>	Date:01/09/14 rev.10
		Issued by: QA/QM
		Page 6 of 10
Family	Single and dual Entry I.V. Catheters	
Device	Invasive, surgical type medical device, non active, for short-term use (more than 60 minutes, less than 30 days) (REF. 93/42/EEC/MDD AND 47/2007 EEC/MDD – APPENDIX IX) — Disposable	
Classification	<b>Class II a (REF. MDD 93/42/CEE – APPENDIX IX) – Rule 7</b> <b>Class II (RIF. CMDR SOR 98-282 – SCHEDULE 1) – Rule 1</b> <i>UMDNS: 10727 – Catheters, Intravenous, Peripheral</i> <i>CND (D.M. 20/02/2007): C0101010102 – Aghi cannula senza valvola d'iniezione</i> <i>C0101010101-Aghi cannula con valvola di iniezione</i> <i>COD. CATEG. REF. prEN1874: 10 – Disposable devices</i>	
Definition	<b>Short-term peripheral venous access.</b> Device for peripheral intravenous catheterisation. Only in connection with other medical devices such as infusion set and/or needleless syringes containing a solution it permits the administration.	
Duration of contact	In order to reduce the risk of thrombophlebitis and bacterial invasion in the CVP <i>in situ</i> in the short term, it is recommended that the device stays in the vein 72-96 hours maximum. (Ref. Atlanta CDC)	

### CHECKS ON FINISHED PRODUCT

By statistical sampling, a series of visual, physical and functional checks are done on the finished product, in accordance with the internal procedures, which are divided up as follows:

- conformity check to the product specification;
- packaging and labelling check;
- visual check of the product;
- check of the distance between the end of the catheter and the beginning of the needle
- penetration test
- Pull test between catheter and catheter hub
- Pull test between needle and needle hub
- Check of connections
- Injection valve leakage test (for dual entry catheters)
- Leakage test catheter hub with pressurised water

### CHEMICAL AND BIOLOGICAL CHECKS

An external laboratory of analysis that has been authorised and approved performs the following chemical and biological checks:

- periodic sterility check: on a batch of packaging;
- biological indicators inactivation check: on every sterilisation batch;
- microbiological contamination check on products before sterilisation (Bioburden Test): periodically;
- checking the absence of abnormal toxicity: done periodically on the product;
- checking the presence of bacterial endotoxin (LAL Test): periodically on manufactured batch;
- chemical checks (reducing substances, acidity or alkalinity, residue on evaporation, spectrophotometric examinations): periodically;
- sterilisation residue check (EO, ethylene chloride, ethylene glycol) in products: done periodically and during the revalidation of the sterilisation process.

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## DISPOSITIF MEDICAL

DELTA MED SpA Viadana (MN)	<b><u>MEDICAL DEVICES — TECHNICAL FORM</u></b>	Date:01/09/14 rev.10
		Issued by: QA/QM
		Page 7 of 10

Family	Single and dual Entry I.V. Catheters
Device	Invasive, surgical type medical device, non active, for short-term use (more than 60 minutes, less than 30 days) (REF. 93/42/EEC/MDD AND 47/2007 EEC/MDD – APPENDIX IX) — Disposable
Classification	Class II a (REF. MDD 93/42/CEE – APPENDIX IX) – Rule 7 Class II (RIF. CMDR SOR 98-282 – SCHEDULE 1) – Rule 1 UMDNS: 10727 – Catheters, Intravenous, Peripheral CND (D.M. 20/02/2007): C0101010102 – Aghi cannula senza valvola d'iniezione C0101010101-Aghi cannula con valvola di iniezione COD. CATEG. REF. prEN1874: 10 – Disposable devices
Definition	<b>Short-term peripheral venous access.</b> Device for peripheral intravenous catheterisation. Only in connection with other medical devices such as infusion set and/or needleless syringes containing a solution it permits the administration.
Duration of contact	In order to reduce the risk of thrombophlebitis and bacterial invasion in the CVP <i>in situ</i> in the short term, it is recommended that the device stays in the vein 72-96 hours maximum. (Ref. Atlanta CDC)

CONSTRUCTION MATERIALS		
COMPONENT	MATERIAL	ABBREV.
Peripheral intravenous and intravascular catheter	Polyurethane	PUR
Rx-opaque line	Pharmaceutical barium sulphate (white) – No. 3 lines.	BaSO <sub>4</sub>
Catheter Hub + Caulking pin	Polypropylene + Polyoxymethylene Copolymer	PP + POM
Injection valve	Silicone rubber	/
Port closure	Polypropylene	PP
Introducer needle	AISI 304 stainless steel	INOX
Needle hub	Styrene-Butadiene Copolymer	ABS
Protector	Low-density polyethylene	PE BD
Vent Fittine Body	Styrene-Butadiene Copolymer	ABS
Vent Fitting Filter	High-density polyethylene	PE AD
Male luer lock cap	Polypropylene	PP

THE DEVICES LISTED IN THE PRESENT TECHNICAL SHEET DON'T CONTAIN
<ul style="list-style-type: none"> <li>Natural latex</li> <li>Phthalates</li> </ul> <p>Furthermore It is not foreseen the use of inks imprinted directly on the parts intended for contact with the patient nor any data dry imprinted.</p>

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## DISPOSITIF MEDICAL

DELTA MED SpA Viadana (MN)	<b><u>MEDICAL DEVICES — TECHNICAL FORM</u></b>	Date:01/09/14 rev.10
		Issued by: QA/QM
		Page 8 of 10
Family	Single and dual Entry I.V. Catheters	
Device	Invasive, surgical type medical device, non active, for short-term use (more than 60 minutes, less than 30 days) (REF. 93/42/EEC/MDD AND 47/2007 EEC/MDD – APPENDIX IX) — Disposable	
Classification	Class II a (REF. MDD 93/42/CEE – APPENDIX IX) – Rule 7 Class II (RIF. CMDR SOR 98-282 – SCHEDULE 1) – Rule 1 UMDNS: 10727 – Catheters, Intravenous, Peripheral CND (D.M. 20/02/2007): C0101010102 – Aghi cannula senza valvola d'iniezione C0101010101-Aghi cannula con valvola di iniezione COD. CATEG. REF. prEN1874: 10 – Disposable devices	
Definition	<b>Short-term peripheral venous access.</b> Device for peripheral intravenous catheterisation. Only in connection with other medical devices such as infusion set and/or needless syringes containing a solution it permits the administration.	
Duration of contact	In order to reduce the risk of thrombophlebitis and bacterial invasion in the CVP <i>in situ</i> in the short term, it is recommended that the device stays in the vein 72-96 hours maximum. (Ref. Atlanta CDC)	

### WARNING

Preventing the infections associated with the use of the intravascular catheters requires that strictly aseptic techniques and manoeuvres are adopted and observed during all the stages required: from the preparation (of the material, patient and personnel), to the insertion of the catheter, and the management of the catheter until it is replaced or removed.

Reconditioning and/or reusing the medical device – although not foreseen- could risk altering the device's performance and be dangerous due to contamination and infection caused by prior use by another patient.

### PACKAGING AND STORAGE CONDITIONS

The primary packaging of the product is a thermoformed and thermo-sealed blister pack made of medical grade paper and plastic film.

The secondary packaging consists of a small cardboard box in cellulose, containing 50 pieces.

The transport packaging consists of corrugated cardboard of suitable dimensions for the content, containing 10 secondary packages (small boxes) equal to 500 pieces.

The opening of the primary packaging (strip) does not require the use of scissors or cutting tools. A non-sealed piece (*Pee!*) actually allows the external packaging to be opened easily.

The opening of the secondary packaging in cardboard also does not require the use of scissors or cutting tools as it is equipped with a tab that allows the easy opening and closure of the lid of the packaging.

As stated on the secondary packaging, the requirement for sterility is guaranteed when the packaging is integral therefore these cannot be used if there is any damaged packaging given that the product cannot be sterilised again.

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## DISPOSITIF MEDICAL

DELTA MED SpA Viadana (MN)	<b><u>MEDICAL DEVICES — TECHNICAL FORM</u></b>	Date:01/09/14 rev.10 Issued by: QA/QM Page 9 of 10
Family	Single and dual Entry I.V. Catheters	
Device	Invasive, surgical type medical device, non active, for short-term use (more than 60 minutes, less than 30 days) (REF. 93/42/EEC/MDD AND 47/2007 EEC/MDD – APPENDIX IX) — Disposable	
Classification	<b>Class II a</b> (REF. MDD 93/42/CEE – APPENDIX IX) – <b>Rule 7</b> <b>Class II</b> (RIF. CMDR SOR 98-282 – SCHEDULE 1) – <b>Rule 1</b> UMDNS: <b>10727</b> – Catheters, Intravenous, Peripheral CND (D.M. 20/02/2007): <b>C0101010102</b> – Aghi cannula senza valvola d'iniezione <b>C0101010101</b> -Aghi cannula con valvola di iniezione COD. CATEG. REF. prEN1874: <b>10</b> – Disposable devices	
Definition	<b>Short-term peripheral venous access.</b> Device for peripheral intravenous catheterisation. Only in connection with other medical devices such as infusion set and/or needleless syringes containing a solution it permits the administration.	
Duration of contact	In order to reduce the risk of thrombophlebitis and bacterial invasion in the CVP <i>in situ</i> in the short term, <b>it is recommended that the device stays in the vein 72-96 hours maximum.</b> (Ref. Atlanta CDC)	

### PACKAGING AND STORAGE CONDITIONS

Packing	Primary	Thermoformed and thermosealed blister pack, with medical grade mesh paper and plastic film.
	Secondary	Cellulose cardboard box, containing 50 pieces
	For Transport	Corrugated cardboard, containing 10 boxes or 500 pieces
Storage conditions	The packages must be stored in dry surroundings that are sheltered from storms and direct and continuous sunlight, in environments free from steam or toxic substances and which are not subjected to frequent changes in temperature. Avoid overloading the transport packaging.	

### DISPOSAL METHODS

Immediately after being removed, the mandrel-needle must be inserted into a rigid container suitable for the disposal of needles and cutting tools.

The medical device, the subject of this Technical Sheet, *must* be disposed of as hospital waste and *cannot* be sent to the incinerator.

The primary packaging (blister pack) *does not* contain PVC.

### STERILISATION

Type	EO (ethylene oxide)
Mixture	EO/CO <sub>2</sub> ratio 10:90 at Bioster EO/CO <sub>2</sub> ratio 15:85 at SterilMilano
Validation	The sterilisation process is validated in accordance with the harmonised European standard EN 550 and the ISO 11135 standard, as well as with F.U.I. and F.E.

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## DISPOSITIF MEDICAL

DELTA MED SpA Viadana (MN)	<b><u>MEDICAL DEVICES — TECHNICAL FORM</u></b>	Date:01/09/14 rev.10 Issued by: QA/QM Page 10 of 10
Family	Single and dual Entry I.V. Catheters	
Device	Invasive, surgical type medical device, non active, for short-term use (more than 60 minutes, less than 30 days) (REF. 93/42/EEC/MDD AND 47/2007 EEC/MDD – APPENDIX IX) — Disposable	
Classification	Class II a (REF. MDD 93/42/CEE – APPENDIX IX) – Rule 7 Class II (RIF. CMDR SOR 98-282 – SCHEDULE 1) – Rule 1 UMDNS: 10727 – Catheters, Intravenous, Peripheral CND (D.M. 20/02/2007): C0101010102 – Aghi cannula senza valvola d'iniezione C0101010101-Aghi cannula con valvola di iniezione COD. CATEG. REF. prEN1874: 10 – Disposable devices	
Definition	Short-term peripheral venous access. Device for peripheral intravenous catheterisation. Only in connection with other medical devices such as infusion set and/or needleless syringes containing a solution it permits the administration.	
Duration of contact	In order to reduce the risk of thrombophlebitis and bacterial invasion in the CVP <i>in situ</i> in the short term, it is recommended that the device stays in the vein 72-96 hours maximum. (Ref. Atlanta CDC)	

### SHELF LIFE

Shelf Life	The devices listed in this technical form have validity of 59 months.
------------	---

### BIOCOMPATIBILITY (ref. EN ISO 10993)

#### Introduction:

Given that this is an invasive surgical type device with the function of peripheral intravenous catheterisation, only in connection with other medical devices such as infusion set and/or needleless syringes containing a solution it permits the administration. The Catheter has biocompatibility requirements, with particular reference to the biocompatibility tests relative to the contact with the material (PUR) of the Catheter itself.

The duration of the contact is between 24 hours and 30 days. (Ref. EN ISO 10993-1).

#### ➤ Catheter:

The material used conforms to the chemical tests of the F.U.I. - Vol. I "Tubular apparatus for blood transfusion and fraction".

The biocompatibility tests performed on the Catheter were done in accordance with the harmonized series of standards EN ISO 10993 and were all passed with positive result.

#### ➤ Needle:

Conforms to the standard ISO 10555-5 and EN ISO 9626 for the pertinent aspects (smooth surfaces, free from roughness, cleaning), being a guide-needle and not a needle that is intended for the transfer of liquids and/or hypodermic needles.

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**Delta Med SpA**