LUXAMED®

LUXAMED® LuxaScope Auris LED

Deutsch	2
English.....	8
Français.....	14
Italiano	20
Español.....	26
Polski	32
Português	38
Suomi	44
EMC Guidelines.....	50

LUXAMED®

LUXAMED GmbH & Co. KG *

Daniel-Weil-Str.3 * 89143 Blaubeuren * Germany

Telefon: +49(0)7344 92905 0 * Fax: +49(0)7344 92905 10

E-Mail: info@luxamed.de * Internet: www.luxamed.de

- DEUTSCH -

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- ENGLISH -

To ensure proper use of this LUXAMED® product, please read this user manual carefully and Retain for future reference

LUXAMED® LuxaScope Auris LED

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1. Otoscope x1
2. Disposable ear tips, grey, each 10 items
2.5 mm
4.0 mm
3. User manual x1
4. Mains adapter for the power supply GTM41078-0605-USB x 1, only for models 3.7 V USB variant
5. Plug: EU / UK / US x 1, only for mains adapter GTM41078-0605-USB, only for models 3.7 V USB variant
6. Lithium-ion battery (10440 350mAh) x1, only for models 3.7 V USB variant
- 7 .Alkaline Batteries (AAA/LR03) x2, only for models 2.5 V variant
8. USB cable x 1, only for models 3.7 V USB variant

Intended Use:

LUXAMED® LuxaScope Auris LED otoscopes are to be used for the examination of the outer ear canal. Furthermore, the instrument can be used for other illumination purposes with non-invasive examinations.

Contraindication:

- Avoid direct contact with mucosa or injured shin.
- Do not shine into the eyes directly.


Warning:

Please make sure that, prior to any use, the ear tips are fixed correctly. NEVER use without ear tips. During examination, ear tip should not remain more than 1 minute in the ear (skin contact).

Inflation for pneumatic test must be done very carefully (risk of excess pressure and any resulting injuries).

-  Use the otoscopes only with LUXAMED® spare parts and accessories.
-  Never use for eye examination.
-  Attention: The slot for fixing the ear funnels is sharp-edged, risk of injury
-  Ensure sufficient contact quality between the otoscope and the charging plug during the charging phase.
-  Do not use the unit while it is charging, do not charge unattended.
-  When cleaning/disinfecting devices that have been used on patients, be sure to use appropriate
-  Personal protective equipment and follow all safety precautions/warnings.
-  Always remove batteries prior to cleaning or disinfecting any device (not rechargeable batteries). Make sure that no liquid penetrates into the housing.
-  Use of any accessories or materials not indicated in the user's manual can degrade the minimum safety of the equipment
-  During use metal components near instrument head can become warm. This is especially true if device is on for extended periods of time. Do not leave device on when not in use.
-  Oscopes are not serviced or maintained while in use with the patient.
-  Regarding any significant RISKS of reciprocal interference posed by the presence of the otoscopes during specific investigations or treatments
-  Regarding potential electromagnetic or other interference between the otoscopes and other devices advice on ways to avoid or minimize such interference.
-  It may happen that the unit switches on unintentionally from the switched-off state due to an electrostatic disturbance. The unit can then be switched off again manually at any time or automatically after 3 minutes.
-  In the diagnostic process, if patients have a burning sensation, user need to stop diagnosing temporarily and turn off the otoscope LED immediately.
-  All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State in which the user and / or the patient is established.

Application:

Insert the mounting socket with detent balls (5) for fixing onto the instrument handle.

Fixing ear tips:

LUXAMED® LuxaScope Auris LED otoscopes are only to be inserted into the ear canal with fixed ear tips (grey). For this, the ear tip can be fixed with the locking pin by slightly turning it clockwise - 12 - through the slot upon the tip of the otoscope (4).

Ear funnels need not be fitted for any other lighting purpose as long as there is no direct contact with the skin or mucous membranes.

The magnifying glass has an approx. 3-fold magnification. The glass can be pivoted in clockwise direction.

Application of pneumatic test:

The LUXAMED® LuxaScope Auris LED otoscope has a connection option for the insufflation bulb for pneumatic test (6).

Power supply:

- 2 Alkali cells (AAA/LR03) or LUXAMED rechargeable battery
- LUXAMED® rechargeable battery with integrated charging electronics. Charging time: Approx. 4 h

The charging lid indicates via the charging status indicator (LED) when the battery is fully charged (LED goes out)
ATTENTION: Fully charge the battery before using it for the first time!

When using the charging cable on a computer it must meet the requirements of IEC 60601-1

Instrument connection:

Connecting port with groove (1).

Switch on:



Press the push button (2) once.

Adjust the brightness:

Adjust the brightness upwards: while the power is on, press and hold the push button (2) until the desired brightness is achieved.

Adjust the brightness downwards:

press and hold the push button (2) again until the desired brightness is achieved.

Switch off:

While the power is on, press the push button (2) once. In addition, this otoscope has an automatic switchoff function after 3 minutes

Changing batteries:

Screw off the battery cap (3) and remove the batteries. Insert the new batteries (AAA/LR03) with the positive pole (+) in direction of the instrument head. Re-screw the battery cap. Never change batteries next to patient.

Rechargeable Battery Replacement:

 The rechargeable batteries may only be changed by authorised personnel. In the event of a defective rechargeable battery, send the handle to LUXAMED® or to an authorised specialist dealer.

Cleaning:

Clean with a soft and lint-free cloth, if necessary, moistened with a little alcohol. The otoscope must not be immersed in liquid.

 Ear funnels must only be used once. Reprocessing may cause damage to the ear.

Disinfection:

A wipe disinfection must be used for disinfection or disinfectant cleaning of the outer surfaces. Disinfectants that are approved for medical devices made of metal, plastic and glass can be used for this purpose. The instructions of the manufacturers of the treatment products must be complied with.

In accordance with the RKI recommendation (Robert Koch Institute, Germany), the disinfectant must fall within the AB range of action (proven bactericidal, fungicidal and virucidal) of the definition of the RKI list

Sterilisation:

Sterilisation is not permitted.



Disposal:

 The product must be disposed of in a separate collection of electrical and electronic equipment.

Technical Specifications:

Nominal voltage	2.5V or 3.7V
Nominal current	23mA
Application part	BF
Protection class	Class II
IP code	IP20
Ambient conditions:	5° C -40° C 10-90% rel. humidity 700-1070 hPa atmospheric pressure
Power supply	2×AAA Alkaline batteries (2.5 V otoscope) 1×3.7V Lithium rechargeable battery (3.7 V otoscope)
Storage / Transport Conditions:	-10° and +45°C. 10% - 90% rel. humidity
Versorgungsspannung	2×AAA Alkaline Batterien (2,5 V Otoskope) 1×3.7V Lithium Akku (3,7 V Otoskope)

Explanation of symbols

	Application part type BF
	Protection class II
	Please observe the Instructions for Use
	Stand-by (On-/off button)
	Do not reuse
	Manufacturer
	Date of manufacture
	Attention
	The product must be disposed of in a separate collection of electrical and electronic equipment.
	Complies with EU directives
	Medical device
	Compliant with UK regulations

Guarantee Conditions

This product is covered by a factory warranty for 2 years from the date of purchase. The date of purchase must be proven by the invoice. During the guarantee period, defects resulting from material or manufacturing faults will be repaired free of charge. The guarantee service shall not extend the guarantee period for the entire device, but only for the replaced components.

Excluded from the warranty is natural wear and tear due to normal use, transport damage and all damage that can be traced back to improper handling (for example, non-observance of the instructions for use) or to tampering with the device by unauthorized persons. The guarantee shall not be grounds for any claims against us for damages. The statutory rights of the purchaser to claims for defects pursuant to § 437 BGB (German Civil Code) shall not be restricted.

In the event of a claim under the guarantee, the device together with the original proof of purchase must be sent to the manufacturer:

LUXAMED GmbH & Co. KG
Daniel-Weil-Straße 3
89143 Blaubeuren
GERMANY

Customer Service:

Guarantee and repair work must be carried out by trained and authorized personnel. The device must not be modified without the permission of the manufacturer. Send the device carefully packed and with properly paid postage to your authorized specialist dealer or directly to:

LUXAMED GmbH & Co. KG
Abteilung Service
Daniel-Weil-Straße 3
89143 Blaubeuren
GERMANY
Tel.: + 49 (0) 7344 92905-0

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EMC Guidelines

Below cables information are provided for EMC reference

Cable	Max. cable length, Shielded/unshielded		Number	Cable classification
Charger USB cable to LuxaScope Auris LED	1m	unshielded	1 Set	DC Power

Important information regarding Electro Magnetic Compatibility (EMC)

 LUXAMED® otoscopes with charging module need special precautions regarding EMC and put into service according to the EMC information provided in the user manual and other documents.

Special precautions need to be observed:

 The LUXAMED® otoscope is only to be used with a medical approved charger, plug and USB cable;

 The use of accessories and cable other than those specified, with the exception of accessories and cables sold by LUXAMED® as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY or decreased LIFESPAN of the charging module.

 Wrong accessories and problems due to EMC can cause an interruption of the essential performance of supplying light for intubation.

 LUXAMED® otoscopes with charging module should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, the LUXAMED® otoscopes be observed to verify normal operation in the configuration in which it will be used.

Guidance and manufacturer' declaration - electromagnetic emissions - for LUXAMED® otoscope with charging module

Emission

Test	Limit	Electromagnetic environment - guidance
Conducted emission	CISPR 11, Group 1, Class B	Device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Radiated emission	CISPR 11, Group 1, Class B	
Harmonic current emissions	IEC 61000-3-2, Class A	Device is suitable for direct connection to the public low-voltage power supply network.

Immunity test levels

Test	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (IEC 61000-4-2)	Contact Discharge: ±8 kV Air Discharge: ±2 kV, ±4 kV, ±8 kV, ±15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF EM field (IEC 61000-4-3)	80-2700 MHz; 1kHz AM 80 %; 10 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>$d = 1.2\sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ for 800 MHz to 2,7 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p>

<p>Proximity fields from RF wireless communications equipment (IEC 61000-4-3)</p>	<p>385 MHz; Pulse Modulation: 18 Hz; 27 V/m 450 MHz, FM + 5 Hz deviation: 1 kHz sine; 28 V/m 710, 745, 780 MHz; Pulse Modulation: 217 Hz; 9 V/m 810, 870, 930 MHz; Pulse Modulation: 18 Hz; 28 V/m 1720, 1845, 1970 MHz; Pulse Modulation: 217 Hz; 28 V/m 2450 MHz; Pulse Modulation: 217 Hz; 28 V/m; 5240, 5500, 5785 MHz; Pulse Modulation: 217 Hz; 9 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance 30 cm.</p>
---	---	--

Electrical fast transients / bursts (IEC 61000-4-4)	Power lines: 2 kV; 100 kHz repetition frequency Signal lines: 1 kV; 100 kHz repetition frequency	Mains power quality should be that of a typical environment.
Surges (IEC 61000-4-5)	L-N: 1kV at 0°, 90°, 180°, 270° L-PE, N-PE: 2 kV at 0°, 90°, 180°, 270°	Mains power quality should be that of a typical environment.
Conducted disturbances induced by RF fields (IEC 61000-4-6)	0.15-80 MHz; 1kHz AM 80 %; 3 Vrms, 6 Vrms in ISM and amateur radio band	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ for 150 kHz to 80MHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).
Rated power frequency magnetic fields (IEC 61000-4-8)	30 A/m, 50 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

<p>Voltage dips / Voltage interruptions (IEC 61000-4-11)</p>	<p>0 % UT for 0.5 cycle at 0° ,45° ,90° ,135° ,180° ,225° ,270° ,315° 0 % UT for 1 cycle at 0° 70 % UT for 25/30 cycles at 0° 0 % UT for 250/300 cycles at 0°</p>	<p>Mains power quality should be that of a typical environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device is powered from an uninterruptible power supply or battery.</p>
--	--	--

www.luxamed.de